

System for quality control of detector liquids during production


Milind Diwan
4/4/2009

Purpose

- Large production facility being planned at Dayabay for various liquids.
- A standardized single system for quality control of optical properties is needed.
- Following proposal is based on the system developed at BNL over last several months.
- BNL system is under use for very precise studies. Johnny Goett has presented results.


Single element

Benefits



- Tube is teflon coated stainless steel. All system components have been certified to be compatible with LAB based LS.
- Optics is based on fibers and LEDs and so it is robust and self calibrating.
- A peristaltic pump (labview controllable) used to obtain samples from production reservoir.
- All components are off-the-self.
- 0.1% stability has been demonstrated.
- A cover gas (for inert atmosphere) can be used.
- Single DAQ can be used for upto 8 tubes.

The entire system



All 5 served by a single VME DAQ system
5 tubes will prevent cross-contamination.
The standard can be used to make sure of
environmental stability.

Capital Cost

	Unit cost	Units	production
SS Tube	\$800	5	3 months
quartz window	\$500	5	2 months
Optics	\$1500	5	1 week
Cables+SMA connectors	\$100	5	1 week
HV	\$1500	3	1 month
VME crate	\$6000	1	1 month
VME controller	\$5000	1	2 months
Flash ADC	\$10000	1	2 months
PMT	\$1500	5	available
Mechanical components	\$5000	1	In house machinist
Total		\$55,500	

We already have PMTs, and the data acquisition system. If needed to save money we can move the BNL system to save ~\$21000

Personnel

- 1 onsite technician for installation 6 months.
- 1 chemistry postdoc
- 1 physics postdoc
- No engineering needed.
- Postdocs needed for both installation and data acquisition.

Space

- 4 meter by 4 meter room is needed with at least 10 ft ceiling. 14 foot ceiling would be very helpful.
- Room needs clean water, storage cabinet for chemicals and, of course clean power.
- Needs to be close to the LS production facility so that samples can be brought in frequently.
- During measurement, which could be for a few hours every night, the room must be isolated from pump vibrations, air conditioning equipment vibrations, and vehicle noise.
- If a section of the room is light-tight, it will make installation easier, but it is not critical.

Summary

- Will require 3-4 months of effort on software.
- Can be assembled quickly.
- Almost all parts are off-the-shelf.
- Will require a technician and postdocs to install.
- No engineering needed.